



Tip Sheet

OMNISSA - EXPORT RESULTS FROM HYPERSPACE


These instructions are for users logging into Hyperspace via the Omnisssa Horizon Client app. End users with Advanced level security (e.g. Quality, Managers, Supervisors, and above) have access to export **Reporting Workbench** and **SlicerDicer Detail Table** results from Hyperspace.

Failure to follow these instructions will result in files being saved in a folder inaccessible to you.


Create A Folder

1. Create a **New Folder** in your OneDrive Documents folder **OR** your Documents folder under This PC.

- a. In your OneDrive – AdventHealth folder




- b. In your Documents folder under This PC




2. Name the new folder **Epic Exports**

Map Folder in Omnissa Horizon Client app


1. Locate the Omnissa icon in your notifications tray at the bottom right of your screen



2. Right click to select **Settings**




3. Select Drive & Folder Sharing




4. Click **Add**
5. Select the **Epic Exports** folder from your OneDrive or Documents folder.
6. Make sure Share your local files is set to **Off**.

Reporting Workbench Export


1. From the Results Viewer activity, click the **Options** button and select **Export Results**.



2. From the Save Results window, give your file a name and click **Browse**.




3. **Critical:** Select the **Epic Exports** folder. **Only the folder mapped in Drive Sharing is available to you later.** The other folders under “This PC” are not accessible to you after export.




4. Click **Save**.
5. **Create** a password for your exported file (this is required).
6. Click **Accept**.

SlicerDicer Export

- From the Detail Table view, right click on a row in the table and select Export All Detailed Rows to Excel or Export All Summary Rows to Excel as needed.



- Create** a password for your exported file (this is required) and click **Accept**.
- Critical:** Select the **Epic Exports** folder. Only the folder mapped in Drive Sharing is available to you later. The other folders under “This PC” are not accessible to you.



- Name** your file.
- Click **Save**.

Date Modified	What Changes	Modified by:	Revised By:
5/14/2025	New version for VMWare only	Amy Gifford	
11/05/2025	Updated for new Omnisia Horizon Client	Heather Dicks	

